



4717 Broadmoor Ave SE, Ste A

Kentwood, MI 49512

P: (616) 949-0080

F: (616) 949-8333

E: jeff@raspberrymed.com

Standard Care Cart



20"W x 22.25"D x 32.6"H

- Compact mobile storage
- 2 - 6" Drawers (12.75" W x 11.5" D x 6" T)
- 3 - 3" Drawers (12.75" W x 11.5" D x 3" T)
- Soft Close drawer slides
- 5" Locking Casters
- Key lock, push button lock available
- Optional side rails for mounting accessories